



Magnesium alloys for biodegradable implants

Laboratory for Metal Physics and Technology
Department of Materials, ETH Zürich

Victor Wessels
Anja Hänzi, Bruno Zberg,
Jörg Löffler, Peter J. Uggowitzer


Coronary stent explanted after 1 month in mini pig

Metamorphosis of the Mg alloy to Ca-P-compound




Micro-CT and LM-Image
after 28 days in mini pig.

Formation of conversion
layer; Ca-P-compound



Biodegradable implants:

- No removal surgery
- Avoidance of long-term adverse reactions (chronic inflammation)
- Early load transfer
- Applications in pediatrics







Witte *et al.*, Biomaterials, 2005



Magnesium:

- Low corrosion resistance
- High biocompatibility
- Favorable mechanical properties
- Young's modulus similar to bone






Cellphone casing

Pressure sensor

Glass: Splats, ribbons (μm)

Bulk metallic glass (cm)




- Specific strength:
Glasses: 337 MPa/g/cm³
 > **Crystalline alloys: 155 MPa/g/cm³**
- High hardness
- High elastic limit, $\epsilon_{el} \sim 2\%$
 Elastic energy density:
 $\frac{1}{2} \sigma_y \epsilon_{el} = \frac{1}{2} E (\epsilon_{el})^2$
- However: Little plasticity

J. F. Löffler, A. A. Kündig, F. H. Dalla Torre, 'Amorphization, rapid solidification, bulk metallic glass processing and properties' (Ch. 17), in CRC Handbook of Materials Processing (Taylor & Francis).

J. F. Löffler, *Int. J. Met. Res.* .2006

Example: Mg–Zn phase diagram




Above 2.4 at.% Zn:
precipitation of intermetallic phases

Metallic glasses:
No such 'phase-diagram limitations'

Song & Atrens, *Adv. Eng. Mater.* **5** (2003) 837.

Oxide layer

Bulk Mg




	Crystalline	Glass
Structure	Multigrain, precipitates	Homogeneous, single phase
Alloying elements	< 5 at.%	< 50 at.%
Corrosion	- Galvanic - Instable oxide layer	- Non galvanic - Stable oxide layer

-Corrosion properties in metallic glasses
can be tailored

Some Mg-based glasses:



(Gu, Shiflet, *et al.*, JMR, 2005)






Here: $\text{Mg}_{(95-x)}\text{Zn}_x\text{Ca}_5$ glasses because of their biocompatibility

Goals

- Produce braided stents from Mg-based glassy wires
- Design Mg-based glasses with minimal degradation-induced H_2 formation for osteosynthesis/fracture fixation




In-house designed melt-extraction setup




→ Production of 100 m long, 100 μm diameter amorphous wires with flawless surface quality, bending plasticity

B. Zberg, E. R. Arata, P. J. Uggowitzer, J. F. Löffler, *Acta Mater.* 2009


Preparation of
dogbone specimens



In-situ tensile testing




Weibull analysis




$$P_i = 1 - \exp \left[- \left(\frac{\sigma_i}{\sigma_{0,Ld^n}} \right)^m \frac{L_i d_i^n}{L_{av} d_{av}^n} \right]$$

Scaling parameter


- Very high tensile strength (4x crystalline alloys) and good reliability
- 4–6% plastic deformation in tension
- Biocompatible MgZnCa wires are interesting for biomedical applications , e.g. self-expandable (braided) stents



Corrosion properties



Oxide layer



Bulk Mg





Mg Glass



	Crystalline	Glass
Structure	Multigrain, precipitates	Homogeneous, single phase
Alloying elements	< 5 at.%	< 50 at.%
Corrosion	- Galvanic - Instable oxide layer	- No galvanic d - Stable oxide layer

Mg-based glasses:

Tailoring of absorption rate
→ Strongly reduced H₂ evolution





Hydrogen evolution during degradation in simulated body fluid



Generally:
 $Mg + 2H_2O \rightarrow Mg(OH)_2 + H_2$



No hydrogen evolution
 above 28 at.% Zn

8 days post-implantation





Correlation of in-vitro and in-vivo hydrogen evolution results

Hydrogen evolution tests




- H₂ evolution decreases by a factor of 100

Linear polarization resistance




- Polarization resistance increases by a factor of 10

- Improved corrosion resistance by Zn addition
- Changing corrosion mechanism due to different surface chemistry



- Impedance $\approx 200 \Omega$
- No change in curve shape with time





- Impedance $> 3,000 \Omega$ after 45 min
- Change in curve shape after 45 min



Buildup of new surface conditions after 45 min


Electrochemical measurements (OCP)

Pourbaix diagram for Zn (calculated)





- Occurrence of elementary Zn
- Zn is in its non-inert state
- Change of surface condition

- Low Zn content:
Slow stabilization of OCP to < -1.0 V (SHE)
- High Zn content:
Fast stabilization of OCP to -0.85 V (SHE)
- Intermediate Zn-content (28 at.%), OCP ≈ -1 V
→ Zn can reach its non-inert state upon pH increase







Hydrogen-evolving
 $Mg_{72}Zn_{23}Ca_5$




- OCP < - 1.0 V:
Zn in its stable metallic form
- Increase of pH:
Zn stays in its stable form

Non-hydrogen-evolving: $Mg_{60}Zn_{35}Ca_5$



- OCP > - 0.9 V:
Zn is always in its non-inert state
→ dense oxide layer

Model:



Implantation of Mg
in various tissues of
domestic pigs



4 tissue types:

- Subcutis
- Muscle
- Liver
- Omentum



2 times:

- 27 days
- 91 days

$Mg_{60}Zn_{35}Ca_5$ glass



Crystal. Mg reference alloy



- All samples show adequate tissue healing reactions, i.e. fibrous capsule foreign-body reaction and neovascularization
- Mg-glasses do not show any *in-vivo* hydrogen evolution

- Mg(95-x)ZnxCa5 metallic glasses were produced and their degradation properties studied in-vitro and in-vivo
- A detailed model of the degradation mechanism was developed based on electrochemical, structural, and spectroscopic analysis of samples for various Zn-content
- Alloys with Zn-content above 28 at.% displayed no clinically observable hydrogen evolution in-vivo
- Alloy wires were used to produce braided stents, and BMG pins are currently being developed for ostesynthesis/fracture fixation

Mg-based metallic glasses could be the next generation of biodegradable implants!

Acknowledgements

Financial support from

- The Swiss Innovation Promotion Agency (CTI-Project 7616.2 LSPP-LS)
- The “BioCompatible Materials and Applications” project (BCMA) initiated by the Austrian Institute of Technology (AIT)
- The Staub/Kaiser foundation (Switzerland) is gratefully acknowledged.